

# 2008 Power Chair Comparo

## Reading this Comparo




Power wheelchairs have been divided — by their manufacturers — into the Rehab or the Consumer category. Within each category, power chairs are listed in alphabetical order by manufacturer.

Manufacturers were asked to report data and specifications for their chairs. When the category was not applicable (such as when no positioning options were available), or when information was not reported or available, "N/A" is listed.

- "RESNA rating" refers to the chair's RESNA dynamic stability rating.
- "Top speed" is the chair's top speed in miles per hour when moving forward.
- "Transit option" indicates whether or not the power chair is equipped for transport in an automotive vehicle while the user is seated in the chair.
- "Transportable" indicates whether or not the power chair disassembles or folds for easy transportability.
- "Target user" describes the manufacturer's intended user, such as "Senior with strength or stamina issues."
- "Power ELRs" = Power elevating legrests.






Data is listed as reported by the manufacturers. This comparo features as many makes and models as space permitted; some entries were edited to meet space requirements. For complete statistics, please consult the manufacturers using the Source List at the end of this comparo.







Manufacturer/Model	RESNA rating	HCPCS code	Drivetrain	User weight capacity	Seat widths/depths	Overall chair width	Overall chair length	Top speed	Turning radius	Seat-to-floor height	Batteries	Transit option?	Tilt/recline?	Positioning options
IMC-Heartway <b>Sahara</b>	10°	K0838	Rear-wheel drive	450 lbs.	18", 20", 22", 18"	25", 46"	6.3 mph	22.5"	20"	12V 50-amp, 2 pcs.	Yes	45°	N/A	
	4°	K0823, K0835, K0848	Center-wheel drive	250 lbs.	16-20", 16-20"	24", 34"	4.8 mph	19"	16.5" starting	2 U1 batteries	Yes, available	55° tilt	5" elevation	
Invacare Corp. <b>TDX Spree</b>	4°	K0848, K0856	Center-wheel drive	165 lbs.	12-16", 12-16"	24", 34"	4.8 mph	19"	14.5" starting	2 U1 batteries	Yes, standard	45° manual tilt, 55° power tilt	5" elevation	
	9°	K0848, K0856, K0861	Front-wheel drive	250 lbs.	17", 19", 21", 16-22"	24.25"-28.5", 43" (34.5" transport length)	5 mph	40"	17.5"	2 Group 34	Yes	45° tilt, 150° recline	8" seat elevation, power elevating legrests	
Permobil <b>C300 Corpus</b>	Tested to 7.5°	K0884	Rear-wheel drive	300 lbs.	17-23", 16-22"	24.25", 36.5"	6 mph	N/A	18.25" with tilt	2 Group 34	Yes	45° tilt	8" adjustable seat height	

Manufacturer/Model	RESNA rating	HCPCS code	Drivetrain	User weight capacity	Seat widths/ depths	Overall chair width	Overall chair length	Top speed	Turning radius	Seat-to-floor height	Batteries	Transit option?	Tilt/recline?	Positioning options
<b>C500 Corpus</b> 	9°	K0884	Front-wheel drive	300 lbs.	17", 19", 21", 16-22"	24.25"-28.5"; 43" (34.5" transport length)	5.5 mph, 7.5 mph optional	44"	19.5"	2 Group 24	Yes	45° tilt, 150° recline	8" seat elevation, power elevating legrests	
	Tested to 7.5°	K0848 without tilt; K0856 with power tilt	Front-wheel drive	165 lbs.	11-16"; 10-18"	24.25"; 34.5"	5 mph	N/A	16.5" with tilt	2 Group 34	Yes	45° power tilt, 45° manual tilt	8" power adjustable seat height on K0856	
<b>Q6000Z</b> 	Incline 10°, Best 4, Worst 4	K0848, K0849, K0856, K0857, K0861, K0868, K0869, K0877, K0884	Mid-wheel 6	300 lbs.	10-24"; 10-24"	25.875"; 36"	6 mph; 8 mph with high-speed motors	22.5"	As low as 16.625"	2 12-volt, deep-cycle NF-22 (Group 3); 2 12-volt, deep-cycle Group 24 (Group 4)	No	55° tilt, 168° recline	Power tilt, recline, tilt/recline, articulating ELRs, articulating foot platform, 7" seat elevation	
	Incline 10°, Best 4, Worst 4	K0848, K0849, K0856, K0857, K0861, K0868, K0869, K0877, K0884	Rear-wheel drive	300 lbs.	10-20"; 10-20"	25.5"; 33"	5 mph (Group 3 standard); 6 mph (Group 4 standard)	23.5"	As low as 16.75"	2 Group 24 or Group 34	No	50° tilt, 168° recline	Power tilt, power recline, power tilt/recline, power articulating elevating legrests	
<b>Rhythm</b> 	9°	K0848, K0850, K0856, K0858, K0861, K0868, K0877, K0884	Mid-wheel drive	300 lbs., 400 lbs. (HD models)	14-24", 15-24" HD models; 12-24"	25.5"; 39"	6.5 mph (Group 3), 5 mph (HD), 8.5 mph (Group 4)	21"	17-20"	Group 24, 22NF	Yes	50° tilt, "Enhanced" recline	11" lift, elevating legrests, articulating legrests	
	9°	K0868, K0884	Rear-wheel drive	400 lbs.	14-24"; 14-22"	25.75" (14-20" seat widths); 27.75" (22" & 24" seat widths); 39"	6 mph (7 or 8.5 mph with Performance packages)	21"	16.8-20.6"	Group 24 or 22NF	Yes	53° tilt, "Enhanced" power recline	Power legrests	

# Consumer

## 2008 PowerChaircompare

Manufacturer/Model	HCPCS Code	Drivetrain	Weight Capacity	Seat Type	Seat Width/Depth	Chair Width/Length	Transportable?	Top Speed	Turning Radius	Batteries	Positioning Options	Target User
<b>Drive Medical</b> Image GT 	K0823	Mid-wheel drive	300 lbs.	Captain's	N/A; N/A	18"; 19"	No	5.5 mph	22.4"	U1	N/A	Users who do not have high-end spinal injuries, but have mobility issues and need the maneuverability of mid-wheel drive
<b>Golden Technologies</b> Compass 	K0823	Center-wheel drive	300 lbs.	Captain's	16-24" in 2" increments; 16-20" in 2" increments	24.5"; 38.5"	No	5 mph	19"	2 U1 or 2 22NF	4" of seat height adjustment	Consumers
<b>Golden-Field</b> Everest CTR 1000 	K0823	Front-wheel drive	300 lbs.	Captain's	19"; 19"	25"; 39"	No	5 mph	20"	2 U1 or Group 22	N/A	Frequent user who requires a highly maneuverable power chair with tight turning capability for indoor/outdoor use
<b>IMC-Heartway</b> Rumba SR-P4R 	K0823	Rear-wheel drive	300 lbs.	Captain's	18-20"; 18"	24"; 41"	No	4.5 mph	20.5"	12V34AH (U1), 2 Pcs.	N/A	Users with degenerative disc disease, degenerative joint disease, neuropathy, edema, etc.
<b>Innovation In Motion</b> X4 Extreme 	N/A	4-wheel drive	400 lbs.	N/A	16", 18", 20", 22", 24"; 16", 18", 20"	27.5"; 46"	No	5.5 mph	53"	Group 4, 75-amp/hour	10" seat elevator, power elevating legrest	Active individuals

Manufacturer/Model	HPCS Code	Drivetrain	Weight Capacity	Seat Type	Seat Width/Depth	Chair Width/Length	Transportable?	Top Speed	Turning Radius	Batteries	Positioning Options	Target User		
Invacare Corp.	Pronto M41		K0823	Center-wheel drive	300 lbs.	Medium-back fold down, semi-recline captain's	16-20"; 16-20"	24", 33.7"	No	5 mph	17.75" front, 19.5" rear	2 U1	No	Geriatric consumers
	Pronto M91		K0848, K0849, K0850, K0856, K0858	Center-wheel drive	300 lbs., 400 lbs.	Semi-recline captain's, semi-recline captain's with solid seat pan, rehab seat (solid base with slingback)	16-24", 16-22"	25.9", 39"	No	6 mph	21.5"	2 22NF	55° power tilt	Active, outdoor user who needs faster speeds and all-terrain performance, or heavy-duty user (300-400 lbs.)
Leisure-Lift	Boss 4.5		K0838 (CC model); K0837 (SS model)	Mid-wheel drive	450 lbs.	Captain's or solid seat	18-28"; 18-23"	25.5", 44"	N/A	5 mph	24.4"	2 Group 22 or Group 24 sealed	50° recline	Bariatric consumers
	Boss 6.75		K0855 (CC model); K0854 (SS model)	Mid-wheel drive	675 lbs.	Captain's or solid seat	24-28"; 21-23"	28.5" -48"	N/A	4.5 mph	27"	2 Group 24 sealed	50° recline	Bariatric consumers
Leisure-Lift	RFP4		K0838 (CC model); K0837 (SS model)	Mid-wheel drive	400 lbs.	Captain's or solid seat	18-24"; 18-21"	25", 40.5"	N/A	5 mph	21.5"	2 Group 22 sealed	50° recline	Bariatric consumers
	Vision Sport		K0822 (model P326-25S); K0823 (model P326-25C)	Mid-wheel drive	300 lbs.	Captain's and solid seat pan	20", 18"	24", 34"	Yes	5 mph	20"	12C U1, 2 pcs.	N/A	N/A

## Manufacturer Source List

### Drive Medical

99 Seaview Blvd.  
Port Washington, NY 11050  
(877) 224-0946  
[www.drivemedical.com](http://www.drivemedical.com)

### Golden Technologies

401 Bridge St.  
Old Forge, PA 18518  
(800) 624-6374  
[www.goldentech.com](http://www.goldentech.com)

### Graham-Field Health Products

2935 Northeast Parkway  
Atlanta, GA 30360  
(800) 347-5678  
[www.grahamfield.com](http://www.grahamfield.com)

### IMC-Heartway

13050 Metro Parkway, Bay 5  
Fort Myers, FL 33966  
(866) 464-9779  
[www.imc-heartway.com](http://www.imc-heartway.com)

### Innovation In Motion

900 Growth Parkway  
Angola, IN 46703  
(800) 327-0681  
[www.mobility-usa.com](http://www.mobility-usa.com)

### Invacare Corp.

One Invacare Way  
Elyria, OH 44036  
(800) 333-6900  
[www.invacare.com](http://www.invacare.com)

### Leisure-Lift

1800 Merriam Lane  
Kansas City, KS 66106  
(800) 255-0285  
[www.pacesaver.com](http://www.pacesaver.com)

### Merits Health Products

730 NE 19th Place  
Cape Coral, FL 33909  
(800) 963-7487  
[www.meritshealth.com](http://www.meritshealth.com)

### Permobil

6961 Eastgate Blvd.  
Lebanon, TN 37090  
(800) 736-0925  
[www.permobil.com](http://www.permobil.com)

### Pride Mobility Products





182 Susquehanna Ave.  
Exeter, PA 18643  
(800) 800-8586  
[www.pridemobility.com](http://www.pridemobility.com)

### Quantum Rehab

182 Susquehanna Ave.  
(866) 800-2002  
Exeter, PA 18643  
[www.quantumrehab.com](http://www.quantumrehab.com)

### Sunrise Mobility

7477 East Dry Creek Parkway  
Longmont, CO 80503  
(800) 333-4000  
[www.sunrisemedical.com](http://www.sunrisemedical.com)

Manufacturer/Model	HCPCS Code	Drivetrain	Weight Capacity	Seat Type	Seat Width/Depth	Chair Width/Length	Transportable?	Top Speed	Turning Radius	Batteries	Positioning Options	Target User
Pride Mobility Products <b>Jazzy 600</b> 	K0822, K0823, K0835, K0836, K0849, K0848, K0856, K0857	Mid-wheel 6	300 lbs.	Captain's or captain's seat with solid seat pan, sling seat, static rehab & power tilt	10-20"; 10-20"	23"; 36.5"	No	5 mph	20"	55AH, NF-22 gel cell or AGM	55° tilt	Active or geriatric users who require tight space maneuverability and indoor/outdoor versatility
Pride Mobility Products <b>Jazzy 614HD</b> 	K0824, K0825	Mid-wheel 6	450 lbs.	Captain's or captain's with solid seat pan, sling, rehab static seating	18-24"; 18-24"	25.125"; 35.5"	No	5 mph	20"	55AH, NF-22 gel cell or AGM	N/A	Bariatric clients
Sunrise Medical <b>Prelude</b> 	K0823	Mid-wheel drive	300 lbs.	Low-back or high-back captain's	18"; 17" (low-back) or 18" (high back)	24"; 33" with-out footplate	Yes	4 mph	21"	U1	N/A	Simple design for consumers using the chair indoors
Sunrise Medical <b>Rhapsody</b> 	K0823 (model CB); K0822 (model BB); K0835 (model SB)	Mid-wheel drive	300 lbs.	Captain's (CB); ASAP (BB); ASAP & SPOT (SB)	16"; 18"; 20" (CB); 16-20" in 1" increments (BB & SB)	23.5"; 32.5" without footrests	No	4 mph	21"	U1 50 ahr AGM	SPOT tilt (50°)	SB: Geriatric users needing pressure relief BB: Adjustable frame seat with flexible base CB: Compact size for more active users