



2007 power chair comparo

Reading This Comparo

Power wheelchairs have been divided — by their manufacturers — into the Rehab or the Consumer category. Within each category, power chairs are listed in alphabetical order by manufacturer.

Manufacturers were asked to report data and specifications for their chairs. When the category was not applicable (such as when no positioning options were available), or when information was not reported or available, “N/A” is listed.




- “RESNA rating” refers to the chair’s RESNA dynamic stability rating.
- “Top speed” is the chair’s top speed in miles per hour when moving forward.
- “Transit option” indicates whether or not the power chair is equipped for transport in an automotive vehicle while the user is seated in the chair.
- “Disassembles/folds” indicates whether or not the power chair disassembles or folds for easy transportability.
- “Target user” describes the manufacturer’s intended user, such as “senior with strength or stamina issues.”
- “Power ELRs” = Power elevating legrests

Data is listed as reported by the manufacturers. This comparo features as many makes and models as space permitted; some entries were edited to meet space requirements. For complete statistics, please consult the manufacturers using the Source List at the end of this comparo.

2007 power chair *compare*

rehab

Manufacturer/Model		RESNA rating	HCPCS code	Drivetrain	User weight capacity	Seat widths/depths	Overall chair width	Overall chair length	Top speed	Turning radius	Seat-to-floor-height	Batteries	Transit option?	Tilt/recline?	Positioning options
Heartway	Sahara RT	8°	K0858	Rear-wheel drive	450 lbs.	18-22"/16-20"	25"	46"	N/A	N/A	N/A	N/A	N/A	N/A	N/A
															
Innovation In Motion	X4 Extreme	N/A	N/A	4-wheel drive	400 lbs.	12-24" & 12-20", both in 2" increments	27.5"	46"	5.5 mph	53"	22"	Group 24, 75-amp/hour	No	35°, 50°, 60° tilt; 135° recline	10" seat elevation, ELRs
															
Innovation In Motion	X5 Extreme	N/A	N/A	Mid-wheel drive	400 lbs.	12-24" & 12-20", both in 2" increments	28"	43"	7.5 mph	21.5"	20.75"	Group 24, 75-amp/hour	Yes	35°, 50°, 60° tilt; 135° recline	10" seat elevation, ELRs
															
nvacare Corp.	3G Torque SP	4	K0868, K0884, K0877	Rear-wheel drive	300 lbs.	12-24"	12-22"	28.5-31.5"	6 mph	32" (Short Base)	17.5"	Two 22NF, two Group 24, three 22NF for vent	Yes	55° tilt, 170° recline	7" seat elevation
															

Manufacturer/Model		RESNA rating	HCPCS code	Drivetrain	User weight capacity	Seat widths/depths	Overall chair width	Overall chair length	Top speed	Turning radius	Seat-to-floor-height	Batteries	Transit option?	Tilt/recline?	Positioning options
Invacare Corp.	TDX SP	4	K0848, K0849, K0850, K0856, K0858, K0861	Center-wheel drive	300 lbs.	12-24"; 12-22"	24-25.5"	35.25"	5.8 mph	22" w/Center Mount footboard	16.5-21"	Two 22NF, two Group 24, three 22NF for vent	Yes	55° tilt, 170° recline	7" seat elevation
															
Leisure-Lift	Boss 6NS SS/CC	8°	K0839, K0838	Mid-wheel drive	600 lbs.	21-32"; 21-24"	28.5"	48"	4.6 mph	27"	15-21.3"	Group 22/27 Gel/SLA	Yes	50° tilt	No
															
Leisure-Lift	Boss 6.75 SS/CC	8°	K0855, K0854	Mid-wheel drive	675 lbs.	24-34"; 21-24"	28.5"	48"	4.5 mph	27"	16-21"	Group 22/24 Gel/SLA	Yes	50°	No
															
Leisure-Lift	RFP4 SS/CC	9°	K0824, K0825	Mid-wheel drive	450 lbs.	18-28"; 18-21"	25"	40.5"	5 mph	21.5"	15.5"	Group 22 Gel/SLA	Yes	50° tilt	No
															
Leisure-Lift	Scout M1 SS/CC	8°	K0823	Front-/rear-wheel drive	300 lbs.	18-24"; 18-21"	23.5"	39"	4.7 mph	20.5"	16"	Group U1 Gel/SLA	Yes	50°	5-6" seat elevation
															
Leisure-Lift	Scout SS/CC	8°	K0848, K0829	Mid-wheel drive	300 lbs.	18-24"; 18-21"	25"	39"	5 mph	21.5"	16"	Group 22 Gel/SLA	Yes	50°	5-6" seat elevation
															

Manufacturer/Model		RESNA rating	HCPCS code	Drivetrain	User weight capacity	Seat widths/depths	Overall chair width	Overall chair length	Top speed	Turning radius	Seat-to-floor-height	Batteries	Transit option?	Tilt/recline?	Positioning options
Permobil	C300 Corpus	9°	K0848, K0856, K0861	Front-wheel drive	250 lbs.	17", 19", 21"; 16-22"	24.25"-28.5"	43" (34.5" transport length)	5 mph	40"	17.5"	Two Group 34	Yes	45° tilt, 150° recline	8" seat elevation, power ELR
															
Permobil	C500 Corpus	9°	K0884	Front-wheel drive	300 lbs.	17", 19", 21"; 16-22"	25.75-28.5"	43" (34.5" transport length)	5.5 mph; 7.5 mph optional	44"	19.5"	Two Group 24	Yes	45° tilt, 150° recline	8" seat elevation, power ELR
															
Quantum Rehab	Q6000Z	10°	K0848, K0849, K0856, K0857, K0861, K0868, K0869, K0877, K0878, K0884	Mid-Wheel 6	300 lbs.	10-24"; 10-24"	25.875"	36"	6 mph (8 mph optional)	22.5"	16.625" minimum	Two 12V, deep cycle NF-22 (Group 3); two 12V, deep cycle Group 24 (Group 4)	No	55°	Power: tilt, recline, tilt/recline, 7" seat elevation, ELRs, articulating foot platform
															
Quantum Rehab	Q600	10°	K0835, K0841, K0848, K0849, K0856, K0857, K0861, K0868, K0869, K0877, K0878, K0884	Mid-Wheel 6	300 lbs.	10-20"; 10-20"	23"	36.5"	5 mph (6 mph optional)	20"	17-17.685" (varies by seating system)	Two 12V, deep cycle NF-22	No	55°	Power: tilt, recline, tilt-recline, ELRs, articulating foot platform
															
Sunrise Mobility	Zippie Z-bop	N/A	K0848	Mid-wheel drive	165 lbs.	10-18"; 13-20"	22.5"	38" w/footplate	5.5 mph	20.5"	15.5-19.5"	U1, 50 Ahr	Yes	N/A	N/A
															
Trac About	IRV2000	N/A	K0014	Tracks	350 lbs.	16-24"; 14-22"	27"	45"	3.5 mph	26" (varies by footrest)	22.5" standard	Two AGM Group 24 or 27	Yes	80° tilt	8" seat elevation
															

consumer

Manufacturer/Model		RESNA rating	HCPCS code	Drivetrain	User weight capacity	Seat widths/depths	Overall chair width	Overall chair length	Top speed	Turning radius	Seat-to-floor-height	Batteries	Transit option?	Tilt/recline?	Positioning options
Drive Medical	Image	N/A	K0823	Mid-wheel drive	300 lbs.	Captain's	19"; 19"	27.5"	39"	No	6.5 mph	22.4"	Two U1/36AH	N/A	Long-term users
															
Graham-Field	E&J Everest CTR 1000	N/A	K0823	Front-wheel drive	300 lbs.	Captain's	19"; 19"	25"	28.75"	No	5 mph	20"	Two U1	N/A	Wants a higher-end chair that offers more comfort and maneuverability
															

Manufacturer/Model	RESNA rating	HCPGS code	Drivetrain	User weight capacity	Seat widths/depths	Overall chair width	Overall chair length	Top speed	Turning radius	Seat-to-floor height	Batteries	Transit option?	Tilt/recline?	Positioning options
Graham-Field 	N/A	K0823	Rear-wheel drive	300 lbs.	Captain's	19"; 19"	25"	28.75"	No	5 mph	30"	Two U1	N/A	Will be in chair for most of the day; wants a higher-end chair with more comfort
Graham-Field 	N/A	K0814	Rear-wheel drive	250 lbs.	Captain's	18"; 16"	22"	33"	Yes	3.5 mph	21"	Two 12V 12AH	N/A	Needs to go a short distance for a short time; needs portability
Graham-Field 	N/A	K0822	Rear-wheel drive	300 lbs.	Sling	16", 18" 20", 16"	24", 25", 27"	41"	Yes	6 mph	33"	Two U1	N/A	Wants a little more than a scooter, needs portability, is on the go all day
Heartway 	N/A	K0825	Rear-wheel drive	400 lbs.	Captain's	22"; 16"	24.8"	42.9"	No	5 mph	23"	Two 12V 50Ah	N/A	Senior or disabled person
Heartway 	N/A	K0823	Rear-wheel drive	300 lbs.	Captain's	18", 20"; 16"	24"	39.5"	No	4.5 mph	20.5	12C 34Ah (U1)/two 50Ah (optional)	N/A	Senior or disabled person
Innovation In Motion 	N/A	N/A	4-wheel-drive	400 lbs.	Captain's	16-24" & 16-20" both in 2" increments	27.5"	46"	No	5.5 mph	53"	Group 24, 75-amp/hour	10" seat elevation, ELRs, tilt, recline	N/A

Manufacturer/Model		RESNA rating	HCPCS code	Drivetrain	User weight capacity	Seat widths/depths	Overall chair width	Overall chair length	Top speed	Turning radius	Seat-to-floor-height	Batteries	Transit option?	Tilt/recline?	Positioning options
Innovation In Motion	X5 Frontier		N/A	Mid-wheel-drive	400 lbs.	Captain's	16-24" & 16-20" both in 2" increments	28"	43"	Yes	7.5 mph	21.5"	Group 24, 75-amp/hour	10" seat elevation, ELRs, tilt, recline	Wants chair with aggressive outdoor capabilities, but also wants to use chair indoors
Invacare Corp.	Pronto M41		N/A	Center-wheel drive	300 lbs.	Semi-recline van seat	16", 17", 20"; 16", 17", 18"	24"	33.7" (footboard up)	N/A	5 mph	17.75" front, 19.5" rear	Two U1	N/A	Needs basic mobility inside/outside; incapable of using manual chair for long periods
Pride Mobility Products	Jazzy Select 6		N/A	Mid-wheel drive	300 lbs.	Captain's	16", 18", 20"; 16", 18", 20"	23.5"	35.625"	No	4 mph	20.5"	Two U1 SLA or AGM	N/A	General user w/o specific rehab needs who qualifies for Group 2 chair

	Manufacturer/Model	RESNA rating	HCPCS code	Drivetrain	User weight capacity	Seat widths/depths	Overall chair width	Overall chair length	Top speed	Turning radius	Seat-to-floor-height	Batteries	Transit option?	Tilt/recline?	Positioning options
Sunrise Mobility	Quickie Prelude 	N/A	K0823	Mid-wheel drive	300 lbs.	Low back, captain's	18", 17", 18"	24"	32.5" w/o footplate	Seat disassembles & folds	4 mph	21"	Two U1	N/A	COPD, heart disease, arthritis or single-amputation client
	Quickie Rhapsody 	N/A	K0823	Mid-wheel drive	300 lbs.	Captain's	16", 18", 20" for both	23.5"	38" w/footplate	Seat comes off & folds down	4 mph	20"	Two U1 or 50 Ahr	N/A	Requires rehab seat or single power option as condition progresses; requires base w/full suspension
Teftec	OmegaTrac 1-6 	N/A	K0848, K0862	Front-wheel drive	300-, 450-lbs.	Sling, solid	14-30"	24.5" and up	48" (longer with options)	N/A	8 mph	30.5"	2 Group 24 gel cell	N/A	15- to 80-year-olds with active lifestyles
Trac About	IRV2000 	N/A	K0014	Tracks	350 lbs.	Captain's	18-24"; 18-24"	27"	49" (38" w/footrest up)	No	3.5 mph; 5.5 mph option	26"	Two AGM Group 24 or 27	Manual reclining back 30°	N/A

manufacturers' source list

Drive Medical
99 Seaview Blvd.
Port Washington, NY 11050
(877) 224-0946
www.drivemedical.com

Graham-Field Health
Products
2935 Northeast Parkway
Atlanta, GA 30360
(800) 347-5678
www.grahamfield.com

Heartway
13050 Metro Parkway, Bay 5
Fort Myers, FL 33966
(866) 464-9779
www.imc-heartway.com

Innovation In Motion
900 Growth Parkway
Angola, IN 46703
(800) 327-0681
www.mobility-usa.com

Invacare Corp.
One Invacare Way
Elyria, OH 44036
(800) 333-6900
www.invacare.com

Leisure-Lift
1800 Merriam Lane
Kansas City, KS 66106
(800) 255-0285
www.pacesaver.com

Permobil
6961 Eastgate Blvd.
Lebanon, TN 37090
(800) 736-0925
www.permobil.com

Pride Mobility Products
182 Susquehanna Ave.
Exeter, PA 18643
(800) 800-8586
www.pridemobility.com

Quantum Rehab
182 Susquehanna Ave.
(866) 800-2002
Exeter, PA 18643
www.quantumrehab.com

Sunrise Mobility
7477 East Dry Creek
Parkway
Longmont, CO 80503
(800) 333-4000
www.sunrisemedical.com

Teftec
12450 Network Blvd.
San Antonio, TX 78249
(210) 477-0330
www.teftec.com

Trac About
1801 SE 9th St.
P.O. Box 502
Newton, KS 67114
(800) 458-8616
www.tracabout.com